

Customized Wheelchair:		Date:			
customer (lastname, firstname):					
phone:		mobile:			
eMail:		country:			
street, No.		ZIP-code:	place:		
animal name:		<table border="1" style="display: inline-table; vertical-align: middle;"> <tr> <td style="text-align: center;">m</td> </tr> <tr> <td style="text-align: center;">w</td> </tr> </table> breed:	m	w	weight: <i>kg</i>
m					
w					
Diagnosis:			age:		

Put the animal in an upright position for taking measures. Ask somebody to help you if your animal is handicapped and not able to stay upright. Please measure the most percisely.

For Measuring circumferences you will need a measuring tape or a piece of cord.
To measure the chest breadth, you can take a big clamp or use a caliper.

Alternatively, you can lean your animal against a wall and hold a board or something over the shoulders. Then just measure the distance between the wall and the board = chest breadth.

1 chest breadth

cm

2 bust girth

cm

3 thigh measurement

right	cm	Left	cm
-------	----	------	----

4

cm

5 hip width:

cm

pick-up

shipping date:

remark:



signature of the submitter

Please note: Your customization order will be manufactured as soon as you submit this filled-out Data-Sheet to PFAFF TIER-ORTHOPAEDIE and completed payment in advance (respectively successful PayPal transaction).
With your signature, you confirm the correctness of your provided measures. If you are not sure, please measure a second time.
DIETER PFAFF TIERORTHOPÄDIE will not assume any liability for measuring errors.